Respiratory Time-Resolved 4D MR Imaging for RT Applications with Acquisition Times below One Minute

> Christopher M Rank¹, Thorsten Heußer¹, Andreas Wetscherek¹, Nami Saito¹, Asja Pfaffenberger¹, and Marc Kachelrieß¹

¹ German Cancer Research Center (DKFZ), Heidelberg, Germany



Introduction

- Radiotherapy of targets in the thorax and upper abdomen is challenging due to respiratory patient motion
- 4D CT is limited by radiation dose, such that multiple breathing cycles cannot be imaged on a daily basis
- 4D MRI has been proposed, however it requires long acquisition times of several minutes
- Simple approach: Gating
 - divide motion cycle into certain gates and reconstruct data from each gate separately
 - trade-off between acquisition time and an appropriate artifact and noise level
- Sophisticated approaches:
 - compressed sensing-based image reconstruction
 - motion-compensated image reconstruction

4D gated MR (40 s)





Aim of Work

- To reduce the acquisition time for 4D respiratory timeresolved MRI below 1 min
- Difficulty: suppress streak artifacts and noise effectively
- Solution: combine compressed sensing-based image reconstruction with motion compensation



4D joint MoCo-HDTV¹ Schematic Overview



[1] Rank, Heußer, Buzan, Wetscherek, Freitag, Dinkel, Kachelrieß. 4D respiratory motion-compensated image reconstruction of free-breathing radial MR data with very high undersampling. *Magn Reson Med*, early view online.



Image Reconstruction Cost Function^{1,2}

Cost function:

$$C = \|X_{pc}Sf - p\|_2^2 + \mu \operatorname{HDTV} f$$
raw data fidelity total variation

- X_{pc} : motion phase-correlated forward transform
 - : coil sensitivity profiles
 - : 4D image volume
- *p* : measured raw data
 - : weight

S

μ

- HDTV : spatial and temporal total variation
- The first term optimizes the raw data fidelity
- The second term improves the image sparsity by optimizing the spatial and temporal total variation
- Both terms are optimized in an alternating manner
- The cost function is optimized for the complete 4D volume including all motion phases

[1] Ritschl, Sawall, Knaup, Hess, Kachelrieß. Iterative 4D cardiac micro-CT image reconstruction using an adaptive spatio-temporal sparsity prior. *Phys. Med. Biol.* 2012.
 [2] Rank, Heußer, Buzan, Wetscherek, Freitag, Dinkel, Kachelrieß. 4D respiratory motion-compensated image reconstruction of free-breathing radial MR data with very high undersampling. *Magn Reson Med*, early view online.



Estimation of Motion Vector Fields Cyclic Deformable Registration¹

Motion estimation only between adjacent phases
 all other motion vector fields given by concatenation





- Incorporate additional knowledge
 - a priori knowledge of quasi periodic breathing pattern
 - non-cyclic motion is penalized
 - error propagation due to concatenation is reduced



Data Acquisition and Processing

- Volunteer measurements of free-breathing thorax and upper abdomen at 1.5 T (Magnetom Aera, Siemens Healthcare)
- 3D-encoded gradient echo sequence with radial stack-of-stars sampling
- Radial sampling in read-out plane, Cartesian sampling in slice direction
- Golden angle (≈ 111.25°) radial spacing



 Data were sorted retrospectively into 20 overlapping respiratory motion phases (10% width)



Results of 4D MR Reconstructions

4D reference gated 3D motion average 6 min 51 s 40 s 4D gated 40 s 4D joint MoCo-HDTV 40 s





Adaptive RT Planning

plan



Courtesy of Vania Batista, University Hospital Heidelberg, Germany, and Nami Saito, German Cancer Research Center, Heidelberg, Germany.



Adaptive RT Planning





Adaptive RT Planning



Courtesy of Vania Batista, University Hospital Heidelberg, Germany, and Nami Saito, German Cancer Research Center, Heidelberg, Germany.



Summary and Outlook

- High quality 4D respiratory time-resolved MRI is possible based on acquisition times below 1 min
- Image quality is comparable to measurements with tenfold acquisition time reconstructed with standard methods
- The reconstructed 4D images and the estimated motion vector fields might be employed for RT applications:
 - reliable target delineation of moving targets
 - patient-specific margin or gating window definition
 - adaptive RT planning

Next steps:

- test of geometric accuracy of MR sequence
- reduction of computation time
- application of algorithm in an ongoing MRI-guided RT study



Thank You!

The 4th International Conference on Image Formation in X-Ray Computed Tomography

> July 18 – July 22, 2016, Bamberg, Germany www.ct-meeting.org



Conference Chair Marc Kachelrieß, German Cancer Research Center (DKFZ), Heidelberg, Germany

This work was supported by the Helmholtz International Graduate School for Cancer Research, Heidelberg, Germany.

This presentation will soon be available at www.dkfz.de/ct.

Parts of the reconstruction software were provided by RayConStruct[®] GmbH, Nürnberg, Germany.

