

# The Rotate-Plus-Shift C-arm Trajectory: Complete Data with Less than 180° Rotation

Ludwig Ritschl<sup>1</sup>, Jan Kuntz<sup>2</sup>, Marc Kachelrieß<sup>2</sup>  
02.12.2014


<sup>1</sup>Ziehm Imaging GmbH, Nürnberg, Germany

<sup>2</sup>German Cancer Research Center, Heidelberg, Germany

## Content

- > Introduction
  - > Theory of the short scan ( $180^\circ$  + fan angle)
  - > Clinical use of mobile 3D C-arms
- > Method
  - > The rotate-plus-shift trajectory
  - > Implementation on a mobile C-arm
  - > Image reconstruction using a filtered backprojection-based approach
- > Results
  - > Preclinical evaluation
- > Discussion

## How Much Rotation Range Do I Need in CT?



We need to measure every point inside the field of measurement in a  $180^\circ$  range. This is sufficient for accurate reconstruction!


- > CT reconstruction bases on 2D Radon transform
- > Measure line integrals over  $180^\circ$  for every point inside the field of measurement (FOM).



First EMI CT scanner in 1972 worked this way!

## Why Requires the Short Scan $180^\circ + \text{Fan Angle } \varphi$ ?





Today's CT scanners and C-arms are constructed in fan-beam geometry.



The reason for  $180^\circ + \varphi$  rotation is, that you get  $180^\circ$  in **every** point inside the field of measurement (FOM).

## Clinical Use of Mobile 3D C-arms

- > Mobile C-arms are clinically used since 1954
- > 3D Imaging on mobile C-arms since 2002.
- > Main 3D applications: Spine **surgery**, trauma **surgery**, ENT **surgery**
- > All applications require a combination of 2D and 3D imaging:  
**2D usability should be preserved on 3D C-arms.**




## What Are The Geometric Requirements for 2D C-Arms?

- **We need space:**
  - A **large C opening** in vertical and horizontal direction is needed for patient positioning and patient access.
  - A **small C shape** reduces interference with other equipment in the OR and improves handling.




## What Should a Good 2D Mobile C-Arm Look Like?

- **We need space:**
  - A **large C-arm opening** in vertical and horizontal direction is needed for patient access
  - A **small C shape** reduces interference with other equipment in the OR



### Large Isocentric C


- Large C
- Enough vertical free space
- Large horizontal opening



### Small Non-Isocentric C

- Small C
- Enough vertical free space
- Large horizontal opening

## What About 3D Capabilities?




**Large Isocentric C**

- Large C
- Enough vertical free space
- Large horizontal opening

**180° +  $\varphi$  of rotation**

**~195°**



**Small Non-Isocentric C**

- Small C
- Enough vertical free space
- Large horizontal opening

**180° -  $\varphi$  of rotation**


**~165°**



Create a way to reach  $180^\circ$   
angular coverage in every  
point of the reconstructed  
FOM with  $180^\circ - \varphi$  of  
rotation!

## What Happens at $180^\circ - \varphi$ of Rotation?


**Idea:** Translate these rays to the other side of the field of measurement!



We already get  $180^\circ$  in this lower part of the field of measurement!

Rotate + Shift Trajectory Using a Mobile C-arm

Scan to Scan in Spine Position



**180° angular coverage  
in the whole field of view!**

## Technical Implementation on a Mobile C-arm

- > Requires motorization of the horizontal and vertical axes holding the C
- > Use adaptive collimators to collimate rays outside the field of view.



## 3D Reconstruction using Filtered Backprojection

- > Use standard FDK type filtered backprojection for reconstruction

$$f(\mathbf{r}) = \int d\alpha \int du \underbrace{w(\alpha, \hat{u}, \hat{v}) \frac{|(\mathbf{o} - \mathbf{s}) \cdot \mathbf{u}^\perp|}{((\mathbf{r} - \mathbf{s}) \cdot \mathbf{u}^\perp)^2} R_F \cos(\beta) \cos(\gamma)}_{\text{Weighting functions}} \underbrace{p(\alpha, u, \hat{v}(\alpha, \mathbf{r}))}_{\text{Raw data}} \underbrace{k(\hat{u}(\mathbf{r}, \alpha) - u)}_{\text{Kernel}}.$$

- > Define dedicated redundancy weights for the rotate + shift scan.



## Results: C-Arm for Preclinical Studies


- > 165° angular rotation
- > 15.8° fan angle
- > 25 kW x-ray generator
- > 30 cm x 30 cm aSi flat detector
  
- > Volume size 16 cm<sup>3</sup>
- > ~ 400 projections
- > ~ 60-150 nGy detector entrance dose per projection (depending on application)



## Results: Preclinical Evaluation

### Cervical spine C1-C7

Visualization:  
C1 – C2 screwing




**Cadaver: 84 years, male,  
osteoporotic**

By the courtesy of Weill Cornell Medical College New York

# Results: Preclinical Evaluation

## Pelvis

### Illosacral screw



Fresh cadaver


By the courtesy of University of Texas, Houston



## Conclusion

- > Proposed new C-arm trajectory that enables for 3D imaging using „classic“ 2D C-arms.
- > Technology is driven by clinical needs.
- > Technical feasibility on a mobile C-arm was demonstrated.
- > First preclinical results show large potential in a wide field of clinical applications.

Thanks to the team!



Christof Fleischmann      Ziehm Imaging GmbH  
Eva-Maria Ilg  
Alexander Jooß  
Wolfgang Kränzel  
Stephan Menger  
Dr. Sven Steckmann  
Christina Tille  
Paul Gabriel

Thank you!

Bryan May  
Ron Villane      Ziehm Imaging Inc.

Prof. Dr. Marc Kachelrieß      DKFZ Heidelberg  
Dr. Jan Kuntz      DKFZ Heidelberg



This work was supported by DFG grant KA 1678/11.